

Kansas Administrative Regulations

Department of Health and Environment

Article 4.—Maternal and Child Health

28-4-504. Methods of specimen collection.

(a) The specimen shall be collected using kits provided by the department.

(b) The form provided with the kit shall be completed before collection of the blood specimen.

(c) The outlined circles on the filter paper portion of the kit shall be saturated with blood in the manner specified on the filter paper.

(d) The specimen shall be delivered by carrier or mailed first-class to the laboratory after the blood has dried and not later than 24 hours from time of collection.

(Authorized by K.S.A. 65-101; implementing K.S.A. 1998 Supp. 65-181; effective, T-87-48, Dec. 19, 1986; effective May 1, 1987; amended April 14, 2000.)

***** Authenticated Kansas Administrative Regulation *****